

Vet-LITE™ 3D

LIGHTWEIGHT EASY-TO-USE SPLINTING MATERIAL

What is Vet-Lite 3D?



Vet-Lite 3D is a lightweight splinting material made from 3-Dimensional polyester fabric and biodegradable X-Lite® resin.



It has been designed to create a strong yet lightweight material - making it ideal for the clinician and the animal.

How to use Vet-Lite 3D



Vet-Lite 3D is available in three different thicknesses, each one suitable for different applications:

3mm:

- Ideal for smaller joints and splints. Excellent for small animals

5mm:

- Use when greater strength is needed; ideal for larger animals (e.g. large dogs, sheep, goats, foals)

7mm:

- An extremely strong and rigid material; ideal for large animals (e.g. horses) where maximum rigidity is required



Features & Benefits:

- Eco-friendly - over 50% less plastic than traditional materials
- Breathable yet strong
- Lightweight for better animal comfort
- Remouldable unlimited times
- Use wet or dry heat to mould
- Non-toxic material without isocyanates
- Ability to spot heat areas
- No layering required
- Safe to use without gloves or apron



Vet-LITE™ 3D

PRODUCT SPECIFICATIONS & ORDERING GUIDE



Quick, Easy and Clean Splint Fabrication:

- Rapid heat activation (from 30 seconds to 2 mins)
- Different heating options - dry heat (hairdryer, heat-gun, oven) or wet heat (splint bath, very hot water (at +65C))
- Quick setting - 2-3 mins working time; 10-20 mins fully cured
- Options for small animals (cats, dogs) to large animals (horses)
- Easily remoulded for adjustment and reuse
- Clean and easy to handle - no gloves required



Perfect for Animal Comfort:

- Extremely lightweight but strong
- Highly breathable and ventilated
- Easy to see any wound discharge
- Safe for users and animals

VET-LITE 3D SHEETS	
9896	Vet-Lite-3D Sheet 600 x 600 x 3mm, Beige
9897	Vet-Lite-3D Sheet 600 x 900 x 3mm, Beige
9898	Vet-Lite-3D Sheet 600 x 600 x 5mm, Beige
9899	Vet-Lite-3D Sheet 600 x 900 x 5mm, Beige
9800	Vet-Lite-3D Sheet 600 x 600 x 7mm, Beige
9801	Vet-Lite-3D Sheet 600 x 900 x 7mm, Beige